

Calx



Known as: Calcium Hydroxide Paste

Description: It is premixed high pH Calcium Hydroxide Water based paste with stabilizer. To be used as a root canal dressing and filling material during endodontic treatment of teeth and as a temporary filling material for pulp capping and pulpotomies.

Packing	Product Code
3 g x 1 Syr + 3 Applicator Tips	9101/003

Intended Use / Indication:

To be used as a root canal dressing and filling material during endodontic treatment of teeth and as a temporary filling material for pulp capping and pulpotomies

Properties & Advantages

- Non setting paste water-based paste.
- Easy handling and easy direct application at the intended site.
- Ideal for temporary filling of root canal and during pulpectomy.
- Excellent radiopacity.
- Premixed, ready to use paste.
- Syringe type for easy delivery of the paste directly into the canal.

Shelf life and storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com