

Calahyd-RC & Calahyd +



Known as: Calcium Hydroxide Paste (With Iodoform)

Description: Calahyd-RC/Calahyd + is a Calcium Hydroxide oil-based paste supplied in a multi dose syringe with disposable dispensing tips. Difference in color of the paste is due to use of different radiopacifiers. Device is used for endodontic treatment of root canals & exhibits bacteriostatic effects and has > 11 pH.

Packing	Product Code
1.5 g x 2 Syr + 4 Applicator tips	9103/15X2
1.5 g x 2 Syr + 4 Applicator tips	9151/15X2

Properties & Advantages:

- Non-setting pastes with oil base.
- Easy handling and easy direct application at the intended site.
- Ideal for infected root canal and pulpectomy in deciduous teeth.
- Premixed, ready to use paste.
- Syringe type for easy delivery of the paste directly into the canal.
- Disposable tips for prevention of cross-contamination.

Intended Use: Temporary root canal filling / dressing materials, root canal sealing materials.

Indications:

- As a temporary Root canal dressing or filling material for weeping root canals.
- As a temporary root canal dressing or filling material during apexification.
- As a root canal filling / sealing material after pulpectomy in deciduous teeth.

Storage & Shelf Life:

- Best before 2 years from the date of manufacturing.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com