

# Mixing Pad



**Description:** Maarc Mixing Pads are made of coated paper.

**Dimensions:** W 60mm x L 60mm & W 75mm x L 146mm.

Packing	Product Code
Mixing Pad Small (1N x 50 Sheets)	8004/001/S
Mixing Pad Big (1N x 50 Sheets)	8004/001/B

## Properties & Advantages:

- High quality mixing pad.
- Easy tear-off disposable coated sheets.
- Liquid impermeable surface.
- Smooth finished surface.
- These mixing pads does not let the dental materials adhere on to their surface.

## Intended Use:

- To be used for kneading or mixing of dental materials.



**Manufactured by:**  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarc dental.com