

Plaster of Paris



Description: Dental plaster or type II gypsum, provides the optimal strength for diagnostic dental models.

Packing	Product Code
1 kg, 3 kg, 5 kg, 20 kg	4201/001, 4201/003, 4201/005, 4201/020.
15 kg bag (Slow Set)	4202/015

Fabrication of diagnostic dental models, mounting of articulators, flasking of trial dentures.

Type II - Specification:

- Water/Powder: 58-60ml / 100gm
- Mixing Time: 15 seconds wet
- Pre-Mix Temperature: Pop and Water 25°C
- Working Time: 3 minutes dry
- Setting Time: 5 to 6 minutes
- Hard Set (Initial): 30 minutes
- Wet Compressive Strength: 150 kg/cm²
- Dry Compressive Strength: 300 kg/cm²
- Setting Expansion (Max): 0.10%

Shelf Life:

- Best before 2 years from the date of manufacturing.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com