

# Pumice Powder



**Description:** Micronized pumice powder to provide a smooth finish to the dentures.

Packing	Product Code
500 g, 5 kg pail, 20 kg bag (Pumice Powder - Fine)	5003/500, 5003/005P, 5003/020
500 g, 5 kg pail, 20 kg bag (Pumice Powder - Medium)	5004/500, 5004/005P, 5004/020
500 g, 5 kg pail, 20 kg bag (Pumice Powder - Coarse)	5005/500, 5005/005P, 5005/020

## Properties & Advantages:

Pumice is available in the form of powder that can be mixed with water to make a paste of desired consistency.

## Indications:

- Polishing tooth surfaces during prophylactic treatment, prior to bleaching and following the removal of orthodontic devices.

## Intended Use:

- For cleaning and polishing the surfaces of teeth.

## Shelf Life:

- Best before 4 years from the date of manufacturing.



Manufactured by:

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com