

# Rootfyx



**Known as:** Bio-Ceramic Root Canal Sealer

**Description:** Rootfyx is a calcium silicate-based root canal sealing and filling material. It is a ready-to-use flowable paste system and is used for sealing /filling the root canal after the complete removal of the infected pulp tissue. Rootfyx is radiopaque and is white in color and can be used in aesthetic regions of teeth.

Packing	Product Code
0.5 g x 4 Syr with 15 Applicator tips	9107/0.5x4

## Properties & Advantages:

- Excellent biocompatibility and good adherence to dentin.
- Good miscibility in water.
- Radiopaque.

## Indications:

- To seal and fill the root canals of teeth.

## Intended Use:

- Sealing and filling of root canals during obturation in the secondary dentition with or without the aid of obturation points.

## Shelf Life:

- Best before 2 years from the date of manufacturing.



**Manufactured by:**

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com