

Sealmax +



Known as: Root Canal Sealer with MTA

Description: Sealmax+ is an MTA based Bioceramic root canal sealing material supplied in the form of powder and resin. The powder and resin are mixed to produce a root canal sealing material of excellent sealing capacity and tissue compatibility. As a result of its close adaptation to the walls and hydrophobic nature, Sealmax + allows tight sealing of root canal during obturation.

Packing	Product Code
3g powder, 5ml liquid and dispensing scoop	8011/001

Properties & Advantages:

- **Filled with Bio-ceramics (MTA):** Contains more than 20% of bio-ceramic material.
- **No moisture required for setting:** Activators in powders makes resin set.
- **Calcium hydroxide release on surface when encounters moisture:** Bio-ceramic fillers on surface when exposed to moisture releases ions and this raises pH of exposed surface, this also helps in inhibiting bacterial growth.
- **Favors bone formation:** Ion release enhances bone formation.
- **Eugenol free:** No long-term post operative pain reported till date.
- Grain free.
- Ideal handling properties.
- Excellent radiopacity.

Indications:

- To seal the root canals of teeth.

Intended Use:

- Sealing of root canals during obturation in the secondary dentition with or without the aid of obturation points (gutta-percha).



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com

Shelf Life:

- Best before 2 years from the date of manufacturing.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarc dental.com