

Calx-LC



Known as: Light Curing Radiopaque Calcium Hydroxide Paste

Description: Calx-LC is a light cure, fluoride-releasing, radiopaque cavity liner and base material specially formulated for use with composites and conventional restorative materials. It chemically bonds to adhesive primers, composites and other resin-based materials, and it micromechanically bonds to dentin.

Packing	Product Code
3g Syr + 4 Applicator Tips	9106/003

Indications:

- Ideal for indirect pulp capping.
- Provides acid protection during the etch technique.
- Serves as a cavity liner under all filling materials.

Properties & Advantages:

- Quick & direct application
- Resistant to acids
- Promotes secondary dentine formation
- Ensures strong adhesion to composites
- Exhibits low solubility in oral fluids

Shelf life and storage:

- Best before 3 years from the date of manufacturing



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com