

Citron (Gel & Solution)



Known as: Root Canal Irrigating Solution

Description: Citron is a slightly thickened, water-soluble Citric Acid gel/solution endodontic cleanser. Intended for cleansing the root canals, it improves surface adherence of root canals to the sealer allowing for ideal seal. Available in two concentrations for ease of use.

Packing	Product Code
Citron 20%: 3g Syringe + 3 Applicator Tips	3022/G/003
Citron 20%: 20ml Syringe + 10 Applicator Tips + 1 Connector	3022/020
Citron 20%: 100ml Bottle Refill	3022/100
Citron 40%: 3g Syringe + 3 Applicator Tips	3023/G/003
Citron 40%: 20ml Syringe + 10 Applicator Tips + 1 Connector	3023/020
Citron 40%: 100ml Bottle Refill	3023/100

Intended Use:

- Citron, is a root canal cleanser, is intended to etch root canal walls just prior to obturation to allow an optimum seal.

Indications:

- Cleanse the root canals during endodontic treatment of teeth.

Properties & Advantages:

- Effectively dissolves Calcium Hydroxide
- Slightly viscous formula enhances lubrication
- Ideal for cleaning and conditioning prepared root canals
- Removes minerals & smear layer from the walls revealing the dentinal tubuli orifices

Shelf life and storage:

- Best before 2 years from the date of manufacturing



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com