

# Gingival Barrier



**Known as:** Liquid Dam

**Description:** Gingival Barrier is designed to isolate soft tissues during bleaching. The methacrylate-based barrier is designed with adequate strength to maintain a barrier yet easily removable from embrasures and undercuts. The contrasting green colour of the Gingival barrier can be clearly identified on the gingiva and teeth for more effective coverage.

Packing	Product Code
3g Syringe + 3 Applicator Tips	8804/003

## Intended Use:

- Gingival barrier is used to isolate soft tissues during the bleaching/whitening of teeth.
- Seal for Rubber Dam.

## Indications:

- To isolate soft tissues during the whitening of teeth or any hazardous process.

## Properties & Advantages:

- Ideal viscosity – Non-Dripping
- Protect gums
- Improve visibility
- Easy to remove
- Fast-setting

## Shelf life and storage:

- Best before 2 years from the date of manufacturing