

# Pattern Resin



**Known as:** Acrylic Resin for Patterns

**Description:** Pattern Resin is a self-curing, general purpose acrylic resin used for making patterns, for fabricating copings, post and core buildups, lingual and palatal bars, implant attachments and telescopic crowns. It can be used for both direct and indirect pattern making.

Packing	Product Code
50g Powder + 2 x 25ml Liquid + 2 Dropper + 2 Measuring Cup + 1 Brush	8997/001

## Intended Use:

- To make the pattern of:
  - Tapered copying for telescopic crowns and other telescopic crowns.
  - Inlay, palatal and lingual bar, connector and clasp.
  - Maryland bridge.
- Temporary splint for soldering.
- In the mouth for core build up pattern.
- Construction of implant super structures.

## Properties & Advantages:

- Excellent Flow & Adaptability
- Minimal Shrinkage
- Quick Setting Time
- Smooth Surface Finish
- High Strength & Stability
- Easy Handling

## Shelf life and storage:

- Best before 2 years from the date of manufacturing



Manufactured by:  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com