

Sealmax Gold



Known as: Resin Based Root Canal Sealer (MTA Filled)

Description: Sealmax Gold is an MTA filled, resin-based root canal sealing material supplied in the form of paste. The base paste and catalyst paste are mixed to produce a root canal sealing material of excellent sealing capacity and tissue compatibility. As a result of its close adaptation to the walls and hydrophobic nature, Sealmax Gold allows tight sealing of root canal during obturation.

Packing	Product Code
20g Base + 12ml Catalyst + 1 Mixing Pad + 2 Nozzle	8031/001

Intended Use:

- Sealing of root canals during obturation in the secondary dentition with or without the aid of obturation points (gutta-percha).

Indications:

- To seal the root canals of teeth

Properties & Advantages:

- Biocompatible
- High Radiopacity
- With Calcium Ions Release
- Easy to mix
- Excellent Flow

Shelf life and storage:

- Best before 2 years from the date of manufacturing



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com