

Silane



Known as: Coupling Agent

Description: Silane is a mono-component solution for application to porcelain following etching with MAARC HF-Gel (hydrofluoric acid).

Packing	Product Code
1.2ml Syringe + 2 Applicator Tips	3025/1.2

Intended Use:

- Priming the surface of porcelain after etching with HF gel (hydrofluoric acid)

Indications:

- Silane is used as a primer during the repair of fractured dental porcelain and porcelain veneers, crowns, or inlays after etching with HF gel
- Silane is used as a primer after etching orthodontic brackets with HF gel before bonding

Properties & Advantages:

- Applicator Tip for easy & precise application

Shelf life and storage:

- Best before 2 years from the date of manufacturing