

# T-Lute NE



**Known as:** Non-Eugenol Temporary Luting Cement

**Description:** T-Lute NE, is a eugenol free self-curing cement to be used for temporary fixation of dental restorations supplied as two paste system. T-Lute NE is especially suitable for use on patients with sensitivity to Eugenol.

Packing	Product Code
30g Base + 13ml Catalyst + 1 Mixing Pad + 1 Nozzle	8991/001

## Intended Use:

- Temporary fixation of provisional restoration.
- Temporary fixation of crowns, bridges, inlays and onlays.
- Fixation of crowns and bridges onto implant abutments.

## Indications:

- Temporary fixation of crowns, bridges & restorations.

## Properties & Advantages:

- Low Film Thickness
- Strong Adhesion
- Easy Removal
- Biocompatible
- Excellent Sealing Properties
- Optimal Fit
- High Retention
- Minimal Cleaning

## Shelf life and storage:

- Best before 2 years from the date of manufacturing



**Manufactured by:**

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com