

Tray Adeziv



Known as: Tray Adhesive

Description: Maarc Tray Adeziv is used to enhance the adhesion of impression materials to stock metal or plastic (solid or perforated) or custom acrylic impression trays. Certain impression materials may get separated from the impression tray upon removal from the mouth leading to distortion errors. Using the Maarc Tray Adeziv helps in securing the impression firmly to the tray upon removal from the mouth.

Packing	Product Code
Tray Adeziv: 15ml Adeziv Thinner: 10ml	3011/001

Intended Use:

- Tray Adeziv is indicated for application to impression trays to enhance the adhesion of impression materials to intraoral impression trays including perforated and non-perforated stock metal and plastic, custom acrylic trays.

Indications:

- To enhance the adhesion of impression materials to intraoral impression trays.

Properties & Advantages:

- Ensures adhesion between the impression material and the impression tray
- Can be used for all kinds of impression trays & materials (made of metal, polypropylene, or methacrylate)
- Firm Adhesion: Ensures strong and reliable bond between alginate and tray
- Thinner provided to reduce wastage & easy cleaning

Shelf life and storage:

- Best before 2 years from the date of manufacturing



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com