

ZOE Impression Paste



Known as: Impression Material

Description: ZOE Impression paste is a Zinc Oxide Eugenol based dental impression material in paste form, easy to mix with excellent flow and dimensional stability.

Packing	Product Code
Zinc Oxide Paste (125g) + Eugenol Paste (85g)	8995/001

Intended Use:

- Useful for edentulous cases and in relining/remodelling
- Techniques with individual trays for uniform support

Indications:

- Impression in edentulous cases & in relining/remodelling of dentures

Properties & Advantages:

- Smooth Consistency & Easy to Mix
- Enough working time to complete border moulding
- Sufficient body to make up for any minor under extensions
- Registers accurate surface details
- Dimensionally Stable
- Cost-effective

Shelf life and storage:

- Best before 3 years from the date of manufacturing



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com