

ZOE Sealer



Known as: Root Canal Obturation Cement

Description: ZOE Sealer is a radiopaque, non-resorbable, Zinc Oxide Eugenol based root canal sealer. Its composition is well tolerated by the tissues and provides adequate sealing of root canals.

Packing	Product Code
15g Powder + 10ml Liquid + 1 Mixing Pad + 1 Dropper + 1 Dispensing Scoop	8021/001

Intended Use:

- Root canal sealing in association with obturation points (gutta-percha).

Indications:

- To seal the Root canals of teeth in association with the gutta-percha point.

Properties & Advantages:

- Root Canal sealer with wide acceptance
- Excellent Handling
- Ensures tight, leak-proof margins
- Non-Resorbable
- Time tested
- Radiopaque

Shelf life and storage:

- Best before 3 years from the date of manufacturing



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com