

# Ziego+



**Known as:** Intermediate Restorative Material

**Description:** Ziego+ (Reinforced Intermediate Filling material) is a type II (ISO 3107) polymer-reinforced zinc oxide-eugenol restorative material designed for intermediate restorations. Ziego+ may also be used as a base under restorative materials and cements that do not contain resin components.

Packing	Product Code
38g Powder + 14ml Liquid + 1 Mixing Pad + 1 Dropper + 1 Dispensing Scoop	8751/001
40g Powder + 15ml Liquid + 1 Mixing Pad + 1 Dropper + 1 Dispensing Scoop	8751/002

## Intended Use:

- Intermediate (temporary) direct restorative material.
- Base material under non-resin-containing restorations.

## Indications:

- Intermediate filling material for restorative treatment of teeth.

## Properties & Advantages:

- High Compressive Strength
- Easy Application
- Easy Removal
- Biocompatible
- Wear Resistance
- Radiopaque

## Shelf life and storage:

- Best before 3 years from the date of manufacturing



**Manufactured by:**  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com