

# SmartSeal RCS



**Known as:** Bioceramic Root Canal Sealing Material

**Description:** SmartSeal RCS is a Calcium Silicate based root canal sealing and filling material. SmartSeal RCS is a ready to use flowable paste system used for sealing/filling the root canal after complete removal of the infected pulp tissue. SmartSeal RCS is radiopaque and white in color and can be used in aesthetic regions of teeth.

Packing	Product Code
1N x 2g Syringe + 6 Applicator Tips	9985/002

**Intended use:** Sealing and filling of root canals during obturation in the secondary dentition with or without the aid of obturation points

**Indications:** To Seal and Fill the Root Canals of teeth

## Properties & Advantages:

- Excellent biocompatibility and good adherence to dentin.
- Good miscibility in water.
- Radiopaque.

## Shelf life and storage:

- Best before 2 years from the date of manufacturing



**Manufactured by:**  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

☎ +919967806543 | ✉ care@shivaproducts.net

🌐 www.maarcidental.com